

Topic 4 – Practical 2

Experimental determination of refractive index

Safety

After prolonged use the ray boxes get hot – avoid contact.

Apparatus and materials

- semicircular Perspex block
- ray box
- power supply
- protractor
- sheet of paper
- ruler
- set square

Introduction


When a ray of light travels through a medium with refractive index n_1 and meets a boundary with another medium of refractive index n_2 , it is partly reflected and partly refracted. If the incident angle is θ_1 and the refracted angle is θ_2 , then:

$$n_1 \sin \theta_1 = n_2 \sin \theta_2 \quad \text{or} \quad \frac{n_1}{n_2} = \frac{\sin \theta_2}{\sin \theta_1}$$

This relationship is known as Snell's law.

Procedure

- 1 Draw a straight line in the middle of the sheet of paper (this will be the line normal to the boundary). Draw another line perpendicular to the first one.



- 2 Using the protractor mark angles of incidence, θ_1 , of 10° , 20° , 30° , 40° and 50° .
- 3 Place the Perspex block on the sheet in such a way that the straight edge is aligned with the second line and the intercept of the two lines is in the middle of the straight edge (point O). Draw the outline of the block so it can easily be placed at the same position.

- ## Questions

-
- Page 2 of 2